

Rotating Mixer



MX-RD-E



MX-RD-Pro



Stainless steel clip
 Cat No.: 18900196, for 1.5ml centrifuge tube
 Cat No.: 18900197, for 15ml centrifuge tube
 Cat No.: 18900198, for 50ml centrifuge tube



Plastic clip
 Cat No.: 18900327, for 1.5ml centrifuge tube
 Cat No.: 18900328, for 15ml centrifuge tube
 Cat No.: 18900329, for 50ml centrifuge tube

Features:

- * Speed is adjustable.
- * LCD display indicates the speed and time for MX-RD-Pro.
- * The plate's angle is adjustable from 0° to 90° for MX-RD-Pro.
- * Material of the plate: Stainless steel with oxidation coating (MX-RD-E), Stainless Steel (MX-RD-Pro).
- * Material of the clip: Plastic (MX-RD-E), Stainless Steel (MX-RD-Pro).
- * Three clip sizes are selectable, suit for the tubes: 1.5ml, 15ml, 50ml.


Accessories for MX-RD-Pro: 18900160 plate or 18900162 plate can be selected as the standard accessory; 18900161 plate is optional.



Cat. No.: 18900160
 Type 1 plate
 60*1.5ml centrifuge tube



Cat. No.: 18900161(Optional)
 Type 2 plate
 16*15ml centrifuge tube



Cat. No.: 18900162
 Type 3 plate
 8*50ml centrifuge tube




Cat. No.: 18900140
 Plate support rods
 4pcs in a group, use with plates

Technical Parameters:

| Model | MX-RD-E | MX-RD-Pro |
|-----------------------|---|---------------------------|
| Speed Range | 0~80rpm | 10~70rpm |
| Inclining Angle Range | / | 0°~90° |
| Timing Range | / | 1~1199min |
| Operation Mode | Continuous | Continuous or Timing |
| Display | / | LCD |
| Motor Type | DC motor | |
| Consumption | 20W | 40W |
| Power Supply | AC110/220V±10%, 50/60Hz | |
| Standard Accessory | 1 pc Plate, 3 bags Plastic clip (12 pcs/ bag) | 1 pc Plate, 1 bag SS clip |
| External Size(W*D*H) | 220*300*310mm | 210*280*300mm |
| Package Size(W*D*H) | 350*450*420mm | 350*400*400mm |
| Gross Weight | 5kg | 7kg |

Vertical Rotating Mixer



MX-RL-E



MX-RL-Pro

Features:

- * Speed is adjustable.
- * LCD display indicates the speed and time for MX-RL-Pro.
- * Material of the tube holder: Stainless steel with oxidation coating (MX-RL-E), Stainless Steel (MX-RL-Pro).
- * Material of the clip: Plastic (MX-RL-E), Stainless Steel (MX-RL-Pro).
- * Three clip sizes are selectable, suit for the tubes: 1.5ml, 15ml, 50ml.



Stainless steel clip
 Cat No.: 18900196, for 1.5ml centrifuge tube
 Cat No.: 18900197, for 15ml centrifuge tube
 Cat No.: 18900198, for 50ml centrifuge tube



Plastic clip
 Cat No.: 18900327, for 1.5ml centrifuge tube
 Cat No.: 18900328, for 15ml centrifuge tube
 Cat No.: 18900329, for 50ml centrifuge tube

Tube Holder for MX-RL-Pro: 18900142 or 18900144 or 18900145 or 18900147 plate can be selected one as the standard accessory; 18900143 and 18900146 plate are optional.



Cat. No.: 18900142
 For 48*1.5ml centrifuge tube
 Held horizontally



Cat. No.: 18900143 (Optional)
 For 24*15ml centrifuge tube
 Held horizontally



Cat. No.: 18900144
 For 24*50ml centrifuge tube
 Held horizontally



Cat. No.: 18900145
 For 32*1.5ml centrifuge tube
 Held vertically



Cat. No.: 18900146 (Optional)
 For 16*15ml centrifuge tube
 Held vertically



Cat. No.: 18900147
 For 16*50ml centrifuge tube
 Held vertically

Technical Parameters:

| Model | MX-RL-E | MX-RL-Pro |
|----------------------|---|--------------------------------|
| Speed Range | 0~80rpm | 10~70rpm |
| Timing Range | / | 1~1199min |
| Operation Mode | Continuous | Continuous or Timing |
| Display | / | LCD |
| Motor Type | DC motor | |
| Consumption | 20W | 40W |
| Power Supply | AC110/220V±10%, 50/60Hz | |
| Standard Accessory | 1 pc Tube holder, 3 bags Plastic clip (12 pcs/ bag) | 1pc Tube holder, 1 bag SS clip |
| External Size(W*D*H) | 530*150*190mm | 510*220*260mm |
| Package Size(W*D*H) | 800*350*420mm | 800*350*420mm |
| Gross Weight | 10kg | 11kg |