

# Corning® Bar Coded Storage Tubes

A system solution for Automated Sample Management


Corning Life Sciences is proud to provide you with the Corning 1D/2D Bar Coded Storage tube line. Compared to the competition, the 1D/2D storage line is a superior storage solution designed to provide maximum identification.

Features include:

- ▶ **MAXIMUM IDENTIFICATION** – Synchronized 2D and linear bar code – along with marking spot
- ▶ **MAXIMUM INFORMATION** – 14 x 14 dot 2D bar code
- ▶ **BAR CODE STABILITY** – Laser-etched, not an attached label
- ▶ **COMPOUND COMPATIBILITY** – Polypropylene construction throughout, making it very inert

## Product Information

	1.3 mL 1D/2D Tube	384 Round Tube	384 Square Tube
Tube material	Polypropylene	Polypropylene	Polypropylene
Temperature range	-180° to +110°C (screw cap) -80° to +110°C (TPE cap)	-80° to +110°C (cap)	-80° to +110°C (cap)
Tube length	43.9 mm	14.0 mm	14.0 mm
Tube height	52.7 mm with screw cap 47.5 mm with TPE cap 53.7 mm total height with tubes and caps, no lid 58.3 mm total height with tubes, caps and lid	16.9 mm with cap 17.2 mm total height with tubes, caps, no lid 23.1 mm total height with tubes, caps and lid	- 17.5 mm total height with tubes, heat seal, no lid 23.1 mm total height with tubes, heat seal and lid
Volume	1.3 mL	120 µL	130 µL
Working volume	1.0 mL	100 µL	110 µL
1D codes	Linear, human readable	–	–
2D code	14 x 14 dots, 14 digits	14 x 14 dots, 14 digits	14 x 14 dots, 14 digits



## Corning® 96 1D/2D Bar Coded Storage Tubes Ordering Information

Cat. No.	Product Description	Qty/Pk	Pk/Cs	Qty/Cs
8500	96 1D/2D Bar Coded Storage Tubes, 1.3 mL, without screw caps, bulk pack	96/pk	10 pk/cs	960 tubes/cs
8501	96 1D/2D Bar Coded Storage Tubes, 1.3 mL, without screw caps, racked with cover	96/rack	10 racks/cs	960 tubes/cs
8502	96 1D/2D Bar Coded Storage Tubes, 1.3 mL, with screw caps, bulk pack	96/pk	10 pk/cs	960 tubes/cs
8503	96 1D/2D Bar Coded Storage Tubes, 1.3 mL, with screw caps, racked with cover	96/rack	10 racks/cs	960 tubes/cs
8504	96 Storage Tube Screw Cap, polypropylene with O-Ring, bulk pack	96/pk	50 pk/cs	4800 caps/cs
8505	96 Storage Tube Screw Cap, polypropylene with O-Ring, on mat, 96 caps per mat, clear with cover	480/pk	10 pk/cs	4800 caps/cs
8506	Thermoplastic elastomer (TPE) 96 Storage Tube Septum Cap, bulk pack	96/pk	100 pk/cs	9600 caps/cs
8507	Thermoplastic elastomer (TPE) 96 Storage Tube Septum Cap on mat, 96 caps per mat	960/pk	5 pk/cs	4800 caps/cs
8508	Empty racks with lids for 96 Screw Capped or Thermoplastic Elastomers (TPE) Capped Storage Tubes	5 racks/pk	2 pk/cs	10 racks/cs
8509	384 2D Bar Coded Storage Tubes, round, without plug caps, racked	384/pk	20 pk/cs	7680 tubes/cs
8510	384 2D Bar Coded Storage Tubes, square, racked	384/pk	20 pk/cs	7680 tubes/cs
8511	384 Storage Tube Plug Cap for round tube, bulk pack	384/pk	100 pk/cs	38400 caps/cs
8512	384 Storage Tube Plug Cap on mats for round tube, 384 caps per mat	3840/pk	10 pk/cs	38400 caps/cs
8513	Empty racks with lids for 384 round or square Storage Tubes	1 rack/pk	20 pk/cs	20 racks/cs

For additional product or technical information, please visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Customers outside the United States, call +1.978.442.2200 or contact your local Corning sales office listed below.

# CORNING

### Corning Incorporated Life Sciences

Tower 2, 4th Floor  
900 Chelmsford St.  
Lowell, MA 01851  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

### Worldwide Support Offices

**ASIA / PACIFIC**  
**Australia/New Zealand**  
t 0402-794-347

**China**  
t 86-21-5467-4666  
f 86-21-5407-5899

**India**  
t 91 124 4604000  
f 91 124 4604099

**Japan**  
t 81 3-3586 1996  
f 81 3-3586 1291

**Korea**  
t 82 2-796-9500  
f 82 2-796-9300

**Singapore**  
t 65 6733-6511  
f 65 6861-2913

**Taiwan**  
t 886 2-2716-0338  
f 886 2-2716-0339

### EUROPE

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 31 20 655 79 28  
f 31 20 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

### All Other European Countries

t 31 (0) 20 659 60 51  
f 31 (0) 20 659 76 73

### LATIN AMERICA

**Brasil**  
t (55-11) 3089-7419  
f (55-11) 3167-0700

**Mexico**  
t (52-81) 8158-8400  
f (52-81) 8313-8589