

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 1395-45DC

Product Description: Three Hole Delivery Cap, GL45

Component Materials:

Cap - PBT
Nut - PEEK (White), ETFE (Brown)
Insert - PTFE
O-ring - Synthetic Rubber

Attributes:

Product Dimensions - 54 mm OD (Cap)

Autoclaving:

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes.

Sterilization:

Shipped Not Sterile

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:

11 SEPT 2014

Rev No: 1