

**Corning Incorporated**  
**Life Sciences**

**Registered**  
**ISO 9001:2008**

## Product Description

**Catalog Number:** 1003P-25

**Product Description:** Corning® Beaker, Low Form, Reusable Plastic, PFA, Not Sterile, 25 mL

**Component Materials:**  
Perfluoroalkoxy-copolymer (PFA)

**Attributes:**

Capacity - 25 mL, Tolerance:  $\pm 1.25$  mL at any mark (5% of full capacity)  
Product Dimensions - 32 mm Body OD x 51 mm High

**Autoclaving:**

This product can withstand autoclave cycles of 121°C and 15 psi for 15 minutes without change to product physical dimensions and volumetric accuracy.

**Heating / Cooling:**

The upper temperature limit for this plastic material is 250°C.

This product IS NOT suitable for Dry Heat Sterilization as it will melt. This product MUST NOT be heated by any direct, surface contacting heat source.

A lower limit usage temperature has not been established for this product.

**Vacuum / Pressure:**

This product is not designed for and is not rated for vacuum or pressure applications.

**Sterilization:**

Shipped Not Sterile

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

**Revision Date:**

01 MAR 2013

Rev No: 1