

## **TECHNICAL SHEET**

PRODUCT NAME: Anti-CD Toxin A mAb (clone TA38)

Code: MT-16TA38

Clone: TA38

Isotype: IgG<sub>1</sub>

Physical State: Liquid

Host Animal (source): Mouse

**Specificity:** TA38 is specific for C. difficile Toxin A detection.

Cross Reactivity: No cross reactivity with Clostridium difficile Toxin B and

Clostridium difficile GDH.

Purification: Anti-CD Toxin A mAb (clone TA38) has been purified in a multi-

step procedure.

≥95% pure by SDS-PAGE.

**Storage Conditions:** Storage at 2-8°C.

**Health & Safety Information:** This product contains sodium azide as preservative. Although the

amount is 0.09%, this substance is poisonous and hazardous and

it should be handled by trained staff only.

Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance

with the respective national regulations.

**Presentation:** PBS pH: 7.4, 0.09% Sodium Azide (NaN<sub>3</sub>).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 14/01/2019



www.certest.es