

# Corning® Step-R™ Repeating Pipettor and Syringe Tips

CORNING

## Corning Step-R

The Corning Step-R is a repeating pipettor suitable to dispense volatile, dense, viscous, or foaming liquids, and for all applications where long dosing series are required. It works with 8 syringe tip sizes to dispense doses from 1  $\mu$ L up to 5 mL. The lightweight construction, ergonomically shaped handle with a finger hook allows for single-handed volume setting, fast syringe tip filling, and convenient serial dispensing. The Step-R has a precise stepping mechanism to ensure repeatability of all dispensed doses, and the positive displacement allows for contamination-free operation.

Each device is tested for functionality and precision according to EN ISO 8655 and is supplied with an individual QC certificate.

## Features

- ▶ Easy volume setting: Volume selection chart conveniently located on the device
- ▶ Fast and precise multiple dispensing: Up to 48 pipetting doses without refilling the tip syringe
- ▶ Wide range of polypropylene syringe tips: 8 different sizes, sterile and non-sterile versions available
- ▶ Positive displacement syringe tips: Dedicated to work with high density or vapor pressure liquids such as glycerol or ethanol
- ▶ Maintenance-free

## Corning Step-R Starter Pack

In order to enable a smooth start, the Corning Step-R Starter Pack is conveniently packaged with a Corning Step-R Repeating Pipettor, 2 syringe tip adapters (25 mL and 50 mL), and 24 syringe tips (0.05 mL, 2.5 mL, 5 mL, 12.5 mL, 25 mL, and 50 mL).



Corning Step-R is supplied with two syringe tip adapters (25 mL and 50 mL).



Corning Step-R Starter Pack

## Ordering Information

### Corning® Step-R™ Repeating Pipettor and Syringe Tips

Cat. No.	Description	Qty/Cs
4217	Corning Step-R Repeating Pipettor	1
6620	Syringe tips, 0.05 mL	100
6628	Syringe tips, 0.05 mL, sterile	100
6621	Syringe tips, 0.5 mL	100
6629	Syringe tips, 0.5 mL, sterile	100
6622	Syringe tips, 1.25 mL	100
6630	Syringe tips, 1.25 mL, sterile	100
6623	Syringe tips, 2.5 mL	100
6631	Syringe tips, 2.5 mL, sterile	100
6624	Syringe tips, 5 mL	100
6632	Syringe tips, 5 mL, sterile	100
6625	Syringe tips, 12.5 mL	100
6633	Syringe tips, 12.5 mL, sterile	100
6626	Syringe tips, 25 mL	25
6634	Syringe tips, 25 mL, sterile	25
6627	Syringe tips, 50 mL	25
6635	Syringe tips, 50 mL, sterile	25

#### Accessories

4218	Adapter for 25 mL syringe tips	1
4219	Adapter for 50 mL syringe tips	1

#### Corning Step-R Starter Pack

Cat. No.	Description	Qty/Cs
6636	Corning Step-R Starter Pack Includes: Corning Step-R Repeating Pipettor (Cat. No. 4217); Adapter for 25 mL syringe tips (Cat. No. 4218) and 50 mL syringe tips (Cat. No. 4219); Syringe tips (5/pk) of 0.05 mL (Cat. No. 6620), 2.5 mL (Cat. No. 6623), 5 mL (Cat. No. 6624), 12.5 mL (Cat. No. 6625); Syringe tips (2/pk) of 25 mL (Cat. No. 6626), and 50 mL (Cat. No. 6627)	1 kit

For more specific information on claims, visit the Certificates page at [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

# CORNING

**Corning Incorporated**  
Life Sciences

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

**ASIA/PACIFIC**  
**Australia/New Zealand**  
t 61 427286832

**China**  
t 86 21 3338 4338  
f 86 21 3338 4300

**India**  
t 91 124 4604000  
f 91 124 4604099

**Japan**  
t 81 3-3586 1996  
f 81 3-3586 1291

**Korea**  
t 82 2-796-9500  
f 82 2-796-9300

**Singapore**  
t 65 6572-9740  
f 65 6861-2913

**Taiwan**  
t 886 2-2716-0338  
f 886 2-2516-7500

**EUROPE**  
CSEurope@corning.com

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 020 655 79 28  
f 020 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

**All Other European Countries**  
t +31 (0) 206 59 60 51  
f +31 (0) 206 59 76 73

**LATIN AMERICA**  
grupoLA@corning.com

**Brasil**  
t 55 (11) 3089-7400

**Mexico**  
t (52-81) 8158-8400