

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 1395-OR

Product Description: PYREX® O-Ring, Silicone Rubber, GL45

Component Materials:
Silicone Rubber, Color: Red

Attributes:
Product Dimensions - 32.7 mm ID x 5.3 mm Thick

Autoclaving:
This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes.

Heating / Cooling:
The normal use temperature range for this product is from -50°C to +125°C. The upper temperature limit for the plastic components of this product is 180°C. This product is not designed for direct contact heating.

Vacuum / Pressure:
This product is a component part and not specifically rated for vacuum or pressure.

Sterilization:
Shipped Not Sterile

TSE/BSE:
Materials of animal origin are not used in the manufacture of this product.

Lot Number Designation:
Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:
30 June 2014

Rev No: 1