

<b>CORNING</b>	<b>Life Sciences</b>	Document Number: <b>LSR00139</b>	Rev.: <b>5</b>
Description (Class): <b>Reference, Customer Technical Data Sheet 60 mm Easy Grip Dish</b>			Page: <b>1 of 2</b>

#### APPLICATIONS:

- Growth of cells.
- Unimpaired observation of dish contents.
- Ideal for microphotography.

#### FEATURES:

- Flat, distortion free optics.
- Lids designed for optimal gas exchange.
- Stacking rings allow for easier stacking and handling.
- Frosted rim improves the handling of the dish. This feature facilitates the ability to pick up the dish without accidentally removing the lid to guarantee aseptic manipulation.
- Reliable vacuum-gas plasma treatment cleans and permanently and consistently modifies the cell-growth surface.
- Non-pyrogenic
- Packaged in peel-open medical-style bags.
- Durable construction for stable dish manipulation

#### REGULATORY COMPLIANCE:

Falcon® products are manufactured under the current ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

#### STERILITY:

Product is gamma irradiated

#### RESTRICTIONS:


- Refer to "Thermoplastics Properties Chart," available at:  
<https://www.corning.com/catalog/cls/documents/compatibility-guides/CLS-DL-AN-226.pdf>

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.

<b>CORNING</b>	<b>Life Sciences</b>	Document Number: <b>LSR00139</b>	Rev.: <b>5</b>
Description (Class): <b>Reference, Customer Technical Data Sheet</b> <b>60 mm Easy Grip Dish</b>			Page: <b>2 of 2</b>

All dimensions in millimeters (mm)

- Assembled Dish and Lid height: 14.94 mm
- Effective growth area = 19.6 cm<sup>2</sup>
- Working Volume = 6.0 - 7.0 mL
- Actual growth area: 19.5 cm<sup>2</sup>



MATERIAL: CRYSTAL-GRADE POLYSTYRENE

ITEM	BDL CAT. NO.	DESCRIPTION	DISH AND LID KEY DIMENSIONS											
			A	B	C	D	E	F	G	H	J	K	L	M
1	353004	60mm Easy Grip Dish. Standard TC. 20/sleeve, 500/case	52.10	50.47	49.86	52.74	53.95	13.13	11.17	.97	54.74	53.11	7.37	.81

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit [www.corning.com/clstrademarks](http://www.corning.com/clstrademarks). All other trademarks are the property of their respective owners.  
© 2013-2019 Corning Incorporated. All rights reserved.