

**Corning Incorporated
Life Sciences**

**Registered
ISO 9001:2008**

Product Description

Catalog Number: 2580

Product Description: Costar® Individual 1 x 8, high binding EIA/RIA strip, flat bottom, without lid

Component Materials:

Strips - Virgin Polystyrene, meets *USP, Class VI* requirements for plastic containers and closures.

Product Dimensions:

Diameter of Well @ top	-	.270 in.	Depth of Well	-	.420 in.
Diameter of Well @ bottom	-	.250 in.	Working Volume per well	-	.36 mL
Tolerances of Dimensions	-	+/- .010 in.			

Surface Characterization:

Surface is characterized to be high binding for EIA/RIA tests. Binding capacity is approximately 500 ng IgG/cm².

Optical Characteristics:

The product is made of clear polystyrene to allow background absorbance to be \pm .005 optical density units from the mean.

Binding Characteristics:

The product has been tested for the attribute of well to well protein binding utilizing a modified ELISA to measure precision. Coefficient of variation must be \leq 8%.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes:	Visual and microscopic examination of the product.
Packaging:	Inspection for accurate labeling, and correct product configuration.
Functional:	Tested for protein binding precision.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 4