

**Corning Incorporated
Life Sciences**

Registered
ISO 9001:2008

Product Description

Catalog Number: 2572

Product Description: Costar ® "egg-crate" flat bottom strip holder

Component Materials:

Holder - Polystyrene, Heavy metal free (meets *CONEG* req.) color concentrate.

Product Dimensions:

Length	-	5.025 in.	Width	-	3.370 in.
Height	-	.560 in.	Tolerances of Dimensions	-	+/- .010 in.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual examination of the product.

Packaging: Inspection for accurate labeling, and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Bovine Spongiform Encephalopathy and Transmissible Spongiform Encephalopathy:

This product complies with the latest revision of EMEA/401/01 "Note for Guidance on minimizing the risk of transmitting animal spongiform encephalopathy agents via human veterinary medicine products" by virtue of all bovine derived material having been processed per specific conditions of section 6.4 of EMEA/401/01.

Rev No: 4