

PRODUCT INFORMATION

Uracil-DNA Glycosylase

Pub. No. MAN0011999

Rev. Date 29 November 2016 (Rev. B.00)

Lot: _ Expiry Date: _

Store at -20 °C

Components	#EN0361	#EN0362	
Uracil-DNA Glycosylase, 1 U/µL	200 U	5 × 200 U	
10X Reaction Buffer	0.4 mL	2 × 1 mL	

www.thermofisher.com

For Research Use Only. Not for use in diagnostic procedures.

Description

Uracil-DNA Glycosylase catalyzes the hydrolysis of the N-glycosylic bond between uracil and sugar, leaving an apyrimidinic site in uracil-containing single or double-stranded DNA. Shows no activity on RNA (1).

Applications

- Control of carry-over contamination in PCR (2).
- Glycosylase mediated single nucleotide polymorphism detection (GMPD) (3).
- Site-directed mutagenesis (4).
- As a probe for protein-DNA interaction studies (5).
- SNP genotyping.
- Cloning of PCR products (6).
- Generation of single strand overhangs of PCR products and cDNA.

Source

E.coli K12 cells.

Molecular Weight

25.6 kDa monomer.

Definition of Activity Unit

One unit of the enzyme catalyzes the release 1 nanomole of uracil from uracil-containing DNA template in 60 min at 37 °C.

Storage Buffer

The enzyme is supplied in: 30 mM Tris-HCl (pH 7.5), 150 mM NaCl, 1 mM EDTA, 1 mM DTT, 0.05% (v/v) Tween 20 and 50% (v/v) glycerol.

10X Reaction Buffer

200 mM Tris-HCl (pH 8.2 at 25 °C), 10 mM EDTA, 100 mM NaCl.

Inhibition and Inactivation

- Inhibitors: Ugi protein from the *Bacillus subtilis* phage PBS2, protein p56 from the *Bacillus subtilis* phage phi29 (7).
- Inactivated by heating at 95 °C for 10 min. Enzyme activity is partially restored at temperatures lower than 55 °C. Therefore put PCR products on ice after PCR and load directly on a gel.

Note

- The abasic sites formed in DNA by Uracil-DNA
 Glycosylase may be subsequently cleaved by heat,
 alkali-treatment or endonucleases that cleave specifically
 at abasic sites.
- UDG is active in the presence or absence of divalent cations.

CERTIFICATE OF ANALYSIS

Endodeoxyribonuclease Assay

No detectable degradation was observed after incubation of supercoiled plasmid DNA with Uracil-DNA Glycosylase.

Ribonuclease Assay

No detectable degradation was observed after incubation of [3H]-RNA with Uracil-DNA Glycosylase.

Labeled Oligonucleotide (LO) Assay

No detectable degradation after incubation of singlestranded or double-stranded radiolabeled oligonucleotides with Uracil-DNA Glycosylase.

Quality authorized by:



Jurgita Zilinskiene

(continued on back page)

References

- 1. Lindahl, T., et al., DNA N-Glycosidases, J. Biol. Chem., 252, 10, 3286-3294, 1977.
- 2. Longo, M.C., et al., Use of uracil DNA glycosylase to control carry-over contamination in polymerase chain reactions, Gene, 93, 125-128, 1990.
- 3. Vaughan, P., McCarthy, T.V., A novel process for mutation detection using uracil DNA glycosylase, Nucleic Acids Res., 26, 810-815, 1998.
- 4. Kunkel, T.A., Rapid and efficient site-specific mutagenesis without phenotypic selection, Proc. Natl. Acad. Sci. USA, 82, 488-492, 1985.
- 5. Devchand, P.R., et al., Uracil-DNA glycosylase as a probe for protein-DNA interactions, Nucleic Acids Res, 21, 3437-3443, 1993.
- 6. Booth, P.M., et al., Assembly and cloning of coding sequences for neurotrophic factors directly from genomic DNA using polymerase chain reaction and uracil DNA glycosylase, Gene, 146, 303-308, 1994.
- 7. Serrano-Heras, G., et al., Protein p56 from the *Bacillus* subtilis phage phi29 inhibits DNA-binding ability of uracil-DNA glycosylase, Nucleic Acids Res, 13, 1-9, 2007.

LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.

The purchase of this product conveys to the buyer the limited, nonexclusive, non-transferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5781 Van Allen Way, Carlsbad, California 92008.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only.* The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Please refer to www.thermofisher.com for Material Safety Data Sheet of the product.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.