

Pressure Infusion Cuffs

Improved "Comfort" with each application!



Pressure Infusion Cuff "Comfort"

Improved "Comfort" with each application!

- For invasive pressure measurement, pressure infusions and pressure transfusions
- Latex-free & suitable for MRI
- Washable cover
- Large tube diameter and large valve outlet for easy and fast pumping up; quick deflating
- Design of the valve-gauge-unit allows a comfortable one-hand operation of bellows and valve
- Optimised fit also for longish infusion bags
- Modular design allows the exchange of single components

Latex-free
&
suitable for MRI

Insert pocket – for recording patient, ward, medication etc.:

- Improves identification and documentation
- Prevents writing on the cover



Pressure Infusion Cuff "Classic"

Quality that pays

- For invasive pressure measurement, pressure infusions and pressure transfusions
- Modular design allows the exchange of single components
- Latex-free versions available
- Washable cover



Pressure Infusion Cuff 500 For single use

- For invasive pressure measurement, pressure infusions and pressure transfusions
- For transfusion/infusion bags up to 500 ml volume
- With pressure gauge
- Cost-efficient and hygienic
- Suitable for single patient use



Pressure Infusion Cuffs • Ordering Information

Comfort-Line – latex free

Order No.	Description	Packaging/Pcs.
97.8100.600	Pressure Infusion Cuff “S-Comfort” - Up to 500 ml volume - With fixed insertion	1
97.8100.610	Pressure Infusion Cuff “K-Comfort” - Up to 1000 ml volume - Variable fit due to velcro	1
97.8100.620	Pressure Infusion Cuff “FR-Comfort” - Up to 1000 ml volume - Solid back and velcro for variable fit - With additional insert pocket for noting down patient, ward, medication etc.	1

Spare Parts

Order No.	Description	Packaging/Pcs.
97.8100.630	Rubber ballon latex-free • Inner tube diameter 8 mm	5
97.8100.631	Valve-Gauge-Unit • Valve-Gauge-Unit for inner tube diameter of 8 mm	1
97.8100.632	Bellows for valve-gauge-unit • Bellows for the connection to the valve-gauge-unit, red, latex-free	1
97.8100.633	Cover for pressure infusion cuff “S-Comfort”	1
97.8100.634	Cover for pressure infusion cuff “K-Comfort”	1
97.8100.635	Cover for pressure infusion cuff “FR-Comfort”	1

Classic-Line

Order No.	Description	Packaging/Pcs.
97.8100.000	Pressure infusion cuff K - Up to 1000 ml volume - Variable fit due to velcro	1
97.8100.011	Pressure infusion cuff K latex-free - Up to 1000 ml volume - Variable fit due to velcro	1
97.8100.100	Pressure infusion cuff Standard - Up to 500 ml volume - With fixed insertion	1
97.8100.111	Pressure infusion cuff Standard latex-free - Up to 500 ml volume - With fixed insertion	1
97.8100.200	Pressure infusion cuff FR - Up to 1000 ml volume - Solid back and velcro for variable fit	1
97.8100.211	Pressure infusion cuff FR latex-free - Up to 1000 ml volume - Solid back and velcro for variable fit	1

Spare Parts

Order No.	Description	Packaging/Pcs.
97.8100.010	Rubber ballon, latex free	5
97.8100.001	Rubber ballon, containing latex	5
97.8100.210	Cover for pressure infusion cuff FR	1
97.8100.002	Cover for pressure infusion cuff K	1
97.8100.102	Cover for pressure infusion cuff Standard	1
97.8100.003	Pressure gauge	5
97.8100.004	Bellows with valve, latex-free	5

Disposable Pressure Infusion Cuff

Order No.	Description	Packaging/Pcs.
97.8100.410	Pressure Infusion Cuff 500 - Up to 500 ml volume - With fixed insertion	2

Distributed by:

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 3050
Fax +49 2293 305-3992
export@sarstedt.com
www.sarstedt.com