

Tuohy needle



- Large, ergonomic, detachable handle plate
- Ideal puncture due to Tuohy bevel
- Centimetre markings for insertion depth control

REGANESTH® Spinal needle

- Reduces post-puncture headache due to atraumatic pencil-point bevel
- Safe handling due to ergonomically designed needle hub
- CSF chamber with magnifying effect enables the identification of even the smallest fluid volumes

Accessories

- Epidural catheter with rounded tip and 3 lateral openings
- Particularly smooth L.O.R. syringe
- Tuohy Borst adapter
- Flat filter 0.2 µm



Ordering Information

Order no.	Description	Pck.
79.8089.159	Epidural set 18G/21Gx900 mm Tuohy needle 18Gx90 mm, epidural catheter 21Gx900 mm, Tuohy Borst adapter, hub insert, flat filter 0.2 µm, L.O.R. syringe LL	5
79.8089.219	CSE set II 18G/27G 900 REGANESTH® Spinal needle pencil-point 27Gx125 mm, Tuohy needle 18Gx90 mm, epidural catheter 21Gx900 mm, Tuohy Borst adapter, hub insert, flat filter 0.2 µm, L.O.R. syringe LL	5
79.8086.000	Flat filter 0.2 µm	20

Distributed by:

SARSTEDT AG & Co. KG
P.O. Box 12 20
D-51582 Nümbrecht
Phone: +49 2293 305 0
Fax: +49 2293 305 3992
export@sarstedt.com
www.sarstedt.com