



PRODUCT INFORMATION

**10X *Taq* Buffer
with KCl**

#B38 4 x 1.25 mL

Lot: _ **Expiry Date: _____**

Store at -20°C

Description

10X *Taq* Buffer with KCl is used in PCR with *Taq* DNA Polymerase.

Does not contain MgCl₂.

10X Buffer Composition:

100 mM Tris-HCl (pH 8.8 at 25°C),

500 mM KCl,

0.8% (v/v) Nonidet P40.

CERTIFICATE OF ANALYSIS

Functionally tested in amplification of 956 bp DNA fragment of a single-copy gene of human genomic DNA.

Quality authorized by:

Jurgita Zilinskiene

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Please refer to www.thermoscientific.com/onebio for Material Safety Data Sheet of the product.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Nonidet is a trademark of Shell Chemical Company. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.