


Revolver Arbeitsrack



- Aufbewahrung unterschiedlicher Probengefäße in einem Rack
- Bohrungen von \varnothing 5 mm - \varnothing 30 mm pro Element
- Elemente sind drehbar und kombinierbar
- Einsetzbar im Temperaturbereich von -80°C bis 121°C
- Abmessungen L x B x H (Element): 108 x 108 x 257 (108 x 108 x 69) mm

Bestellinformation

Bestell Nummer	Bezeichnung	Verpackung
93.879	Revolver Arbeitsrack	1/Beutel 4/Karton

SARSTEDT AG & Co.
Postfach 12 20 · D-51582 Nümbrecht
Telefon +49 22 93 305 - 0
Telefax+49 22 93 305 - 3450
Service 0800 (Deutschland)
Telefon (0800) 0 83 305 - 0
info@sarstedt.com
www.sarstedt.com

Revolver Workstation



Technical modifications reserved

This document may contain information on products that may not be available in particular countries

- Designed to accommodate a range of tube sizes in one rack
- Each unit with Ø 5 mm to Ø 30 mm orifices
- Units are rotatable and freely exchangeable
- Ideal for use in the temperature range from -80°C to 121°C
- L x W x H (per unit): 108 x 108 x 257 (108 x 108 x 69) mm

Ordering Information

Art. No.	Description	Packaging
93.879	Revolver Workstation	1/bag 4/case

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305-0
Fax +49 2293 305-3992
export@sarstedt.com
www.sarstedt.com