


Ordering information

1250 μ l PIPETTE TIPS



LONG & VERSATILE

Advanced

Highly transparent material, clear filling level rings, optimal suitability for multichannel pipettes, the perfect fit – all this is combined in our newly developed extra long tip. It convinces with its precision and comfortable handling together with multi-channel and single-channel pipettes.

Deep-reaching

With a length of 90 mm, our new tip is particularly convenient. Pipetting to the bottom of longer vessels without touching the vessel wall is now even easier. Thanks to the slim design, transfer your sample directly in situ with ease.

Versatile

The stable tray – the perfect connection between pipette and tip and the optimal design of the box. Our newly designed 1250 µl-tip is your reliable tool in the laboratory. Now also available in the low retention version and in all packaging forms.



TOP TIP Discover which is the best SARSTEDT tip for your pipette in our easy and convenient online compatibility table at www.sarstedt.com.

Ordering information

Vol.	Order no.	Description	Filter	Purity	Design	Packaging	Pieces/ carton
1,250 µl	70.3060.275	Low retention filter tip, 1000 µl*, Biosphere® plus, transparent, 96 pcs/Box	+		Low Retention		1,920
	70.3060.375	Low Retention Filter Tip, 1000 µl*, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Low Retention		1,920
	70.3060.255	Filter tip, 1000 µl*, Biosphere® plus, transparent, 96 pcs/Box	+		Standard		1,920
	70.3060.355	Filter tip, 1000 µl*, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Standard		1,920
	70.3060.205	Pipette tip, 1250 µl, Biosphere® plus, transparent, 96 pcs/Box	-		Standard		1,920
	70.3060.305	Pipette tip, 1250 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	-		Standard		1,920
	70.3060.210	Low retention pipette tip, 1250 µl, PCR Performance Tested, transparent, 96 pcs/Box	-		Low Retention		1,920
	70.3060.110	Low retention pipette tip, 1250 µl, PCR Performance Tested, transparent, 5x 96 pcs/StackPack	-		Low Retention		1,920
	70.3060.200	Pipette tip, 1250 µl, PCR Performance Tested, transparent, 96 pcs/Box	-		Standard		1,920
	70.3060.100	Pipette tip, 1250 µl, PCR Performance Tested, transparent, 5x 96 pcs/StackPack	-		Standard		1,920
70.3060	Pipette tip, 1250 µl, PCR Performance Tested, 500 pcs/Bag	-		Standard		5,000	

*Reduced pipetting volume

We print on FSC certified paper.



Further information can be found at www.refillrevolution.tips

SARSTEDT AG & Co. KG
P.O. Box 12 20
D-51582 Nümbrecht
Phone: +49 2293 305 0
Fax: +49 2293 305 3992
export@sarstedt.com
www.sarstedt.com

SARSTEDT