

PRODUCT INFORMATION

25 mM MgCl₂

#AB-0359

5 mL

Lot _

Expiry Date _

Store at -20°C

∐ 74

www.thermoscientific.com/onebio

Description

Magnesium Chloride at 25 mM.

Use

Magnesium chloride is usually used in PCR reactions. Concentration should be optimized according to template and primer combination.

Tips before use

Thaw the reagents on ice. Mix and spin down the solutions before use to recover the maximum amount.

Storage Conditions

Store at -20°C until ready for use. Avoid repeated freeze thawing.

Technical support

EMEA: ts.molbio.eu@thermofisher.com

Americas & APAC: ts.molbio@thermofisher.com

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes* and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.thermoscientific.com/onebio for Material Safety Data Sheet of the product.

© 2012 Thermo Fisher Scientific, Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries