# dNTP/dUTP Mix, molecular biology grade

## Catalog Number R0251

**Pub. No.** MAN0013135 **Rev.** B.00



**WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Sheets (SDSs) are available from **thermofisher.com/support**.

## Contents and storage

Contents	Amount	Storage
dNTP/dUTP Mix, molecular biology grade	1 mL	-25 °C to -15 °C

#### **Description**

dNTP/dUTP Mix is aqueous solution containing dATP, dCTP, dGTP, each at a final concentration of 2 mM and dUTP of 4 mM.

#### **Applications**

For use in PCR, real-time PCR and RT-PCR with *Tag* DNA Polymerase, cDNA synthesis and other applications.

#### Note

dNTP/dUTP Mix is not suitable for PCR with proofreading enzymes or enzyme mixes containing proofreading enzymes.

#### Limited product warranty

Life Technologies Corporation and/or it affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 Vilnius, Lithuania For descriptions of symbols on product labels or product documents, go to **thermofisher.com/symbols-definition** 

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For Research Use Only. Not for use in diagnostic procedures.

 $thermofisher.com/support \mid thermofisher.com/askaquestion$ 

thermofisher.com

