

PrepMan™ Ultra Sample Preparation Reagent

Catalog Number 4318930

Pub. No. 4380077 Rev. C

Note: For safety and biohazard guidelines, see the “Safety” appendix in the *PrepMan™ Ultra Sample Preparation Reagent Protocol* (Pub. No. 4367554). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Additional information and procedures are available in the *PrepMan™ Ultra Sample Preparation Reagent Protocol* (Pub. No. 4367554).

Product information

PrepMan™ Ultra Sample Preparation Reagent provides a simple way to prepare DNA from a wide range of sample types including:

- Processed foods and their ingredients
- Bacteria
- Fungi
- Mammalian tissue smears
- Hair
- Human cells (buccal swab)
- Whole blood
- Plasmid preparations

Contents

Contents	Amount	Storage
PrepMan™ Ultra Sample Preparation Reagent	One bottle, 20 mL	15–25°C on receipt IMPORTANT! Do not freeze or autoclave.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

thermofisher.com/support | thermofisher.com/askaquestion



Manufacturer: Thermo Fisher Scientific | 7 Kingsland Grange | Warrington, Cheshire WA1 4SR | United Kingdom

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Revision history: Pub. No. 4380077 (English)

Revision	Date	Description
C	17 April 2018	Updating template, legal, and contact information.
B	21 January 2014	Baseline for this revision history.

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity: Life Technologies Corporation | Carlsbad, CA 92008 USA | Toll Free in USA 1 800 955 6288

©2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For Research Use Only. Not for use in diagnostic procedures.