

PRODUCT INFORMATION Thermo Scientific 10X *Taq* Buffer without Detergent

#B55 4 x 1.25 mL

Lot: Expiry Date:

Description

10X *Taq* Buffer without Detergent is used in PCR with *Taq* DNA Polymerase. The buffer does not contain detergents which can interfere with certain applications, e.g., microarray experiments. Does not contain MgCl₂.

10X Buffer Composition

100 mM Tris-HCI (pH 8.8 at 25°C), 500 mM KCI.

CERTIFICATE OF ANALYSIS

Tested for the absence of endo-, exodeoxyribonucleases, ribonucleases and functionally tested in amplification of 956 bp DNA fragment of a single-copy gene of human genomic DNA.

Quality authorized by:

Et Jurgita Zilinskiene

Store at -20°C

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only.* The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to <u>www.thermoscientific.com/onebio</u> for Material Safety Data Sheet of the product.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.