

Syringe Labelling

Adhesive labels according to ISO 26825 and the DIM and DGAI recommendations



Syringe labels and label dispensers

- Reduction of substance mix-ups due to standardised colour-coded marking of the drug classes
- Suitable metal label dispenser in two sizes available
- Syringe and perfusor labels

Label database

An overview of the currently available labels can be found using the following link:

<https://www.sarstedt.com/en/products/clinic/syringe-labelling/labels/>

Label dispensers



- Clearly arranged storage
- Quick access
- Easy cleaning
- Suitable as table model and for wall mounting

Colour codes of the drug classes

Hypnotics	Muscle relaxants	Muscle relaxant antagonists
Opiates	Opioid antagonists	Local anaesthetics
Vasopressors	Antihypertensive agents/vasodilators	Benzodiazepines
Benzodiazepine-antagonists	Anticholinergics	Cholinergics
Antiemetics	Different compounds	Anticoagulants
Coagulants	Antiarrhythmics	Anticonvulsants
Bronchodilators	Inodilators	Hormones
Insulin Hormones	Electrolytes	Electrolytes potassium
Electrolytes NaCl		

We comply with the ISO 26825 standard and the DIVI and DGAI recommendations for the manufacture of our syringe labels. This allows us to contribute to a reduction in substance mix-ups and increase patient safety.

The following label variants are available:

Label variants		
Syringe labels Size: 35 x 13 mm, Reel with 1000 pcs.	Perfusor labels Size: 80 x 40 mm, Reel with 500 pcs.	Customized labels Both sizes available Minimum purchase: 6 reels

Ordering information

Order no.	Description
97.8300.500	Label dispenser with cover For 7 reels of syringe labels or 3 reels of perfusor labels
97.8300.501	Label dispenser Maxi with cover For 10 reels of syringe labels or 4 reels of perfusor labels

Distributed by:

SARSTEDT AG & Co. KG P.O. Box 12 20 D-51582 Nümbrecht Phone +49 2293 305 - 0 Fax +49 2293 305 - 3992 export@sarstedt.com www.sarstedt.com
