

# REGANESTH® unoplex NRFit


Ultrasound-guided peripheral nerve blocks  
with ISO 80369-6 NRFit™ Non-Luer Connector

NRFit™




Suitable for ultrasound

*In plane*



*Out of plane*




## REGANESTH® unoplex NRFit

- The REGANESTH® unoplex needle was developed for ultrasound-guided nerve blocks
- The needle tip enables outstanding visibility, even for steep puncture angles
- Reduces the risk of accidental misconnections due to the ISO 80369-6 NRFit™ Non-Luer Connector

## Efficient

The *rhombus cut* ultrasound marker enables excellent echogenicity



## Easy to use

Convenient handling thanks to the optimised and ergonomically designed REGANESTH® handle and bevel position marking



## Accessories



## Ordering information

Order no.	Description	Packaging unit
79.8096.044	<b>REGANESTH® unoplex NRFit 25G x 40 mm</b> with detachable 60 cm injection tube	10
79.8096.054	<b>REGANESTH® unoplex NRFit 22G x 50 mm</b> with detachable 60 cm injection tube	10
79.8096.084	<b>REGANESTH® unoplex NRFit 22G x 80 mm</b> with detachable 60 cm injection tube	10
79.8096.104	<b>REGANESTH® unoplex NRFit 20G x 100 mm</b> with detachable 60 cm injection tube	10
79.8096.254*	<b>Unoplex-Catheter-Set 18G/50 mm NRFit</b> Plexus needle REGANESTH® unoplex 18G x 50 mm, Plexus catheter 21G x 500 mm with stylet, central opening, extension tube 60 cm, Tuohy-Borst adapter, hub insert, flat filter 0.2 µm with fixation, with NRFit™ Non-Luer connector according to ISO 80369-6	5
79.8096.284*	<b>Unoplex-Catheter-Set 18G/80 mm NRFit</b> Plexus needle REGANESTH® unoplex 18G x 80 mm, Plexus catheter 21G x 500 mm with stylet, central opening, extension tube 60 cm, Tuohy-Borst adapter, hub insert, flat filter 0.2 µm with fixation, with NRFit™ Non-Luer connector according to ISO 80369-6	5

\*Available soon

### Distributed by:

SARSTEDT AG & Co. KG  
P.O. Box 12 20  
D-51582 Nümbrecht  
Phone: +49 2293 305 0  
Fax: +49 2293 305 3992  
export@sarstedt.com  
www.sarstedt.com