# Adhesive PCR Foil

Catalog Number AB-0626

Pub. No. MAN0019424 Rev. B.0

## **Product description**

The Adhesive PCR Foil is a strong foil that provides an excellent bond that eliminates evaporation at high temperatures. This makes the foil suitable for both storage purposes (ambient to  $-40^{\circ}$ C) and for PCR. The foil can be easily pierced with a pipette tip, leaving no adhesive build up on the tip.

## Contents and storage

Product	Cat. No.	Amount	Storage <sup>[1]</sup>
Adhesive PCR Foil	AB-0626	100 sheets	<ul><li>Keep cool and dry</li><li>Avoid direct exposure to sunlight</li></ul>

<sup>[1]</sup> It is recommended to use the seals within one year from date of purchase.

## **Product specifications**

	Complete: 130 x 80 mm With the application strip removed: 122 x 80 mm	
Seal thickness	75 µm (without release liner)	
Temperature range	-40°C to 120°C	
	Bond strength is dependent on application pressure, so sufficient pressure must be applied. For optimal seal application, we recommend the use of the Adhesive Seal Applicator (Cat. No. AB-1391).	

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at <a href="http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html">www.thermofisher.com/us/en/home/global/terms-and-conditions.html</a>. If you have any questions, please contact Life Technologies at <a href="http://www.thermofisher.com/support">www.thermofisher.com/us/en/home/global/terms-and-conditions.html</a>. If you have any questions, please contact Life Technologies at <a href="http://www.thermofisher.com/support">www.thermofisher.com/us/en/home/global/terms-and-conditions.html</a>. If you have any questions, please contact Life Technologies at <a href="http://www.thermofisher.com/support">www.thermofisher.com/support</a>.

#### thermofisher.com/support | thermofisher.com/askaquestion



Life Technologies Corporation | 5781 Van Allen Way | Carlsbad, CA 92008

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Important Licensing Information**: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

