



PRODUCT INFORMATION

Thermo Scientific
10X DreamTaq
Buffer

#B65 4 x 1.25 mL

Lot: **Expiry Date:**

Store at -20°C

www.thermoscientific.com/onebio

Description

Thermo Scientific™ 10X DreamTaq™ Buffer is a proprietary buffer optimized for PCR using DreamTaq™ DNA Polymerase. The buffer contains both KCl and (NH₄)₂SO₄ to provide high specificity of primer annealing at a wide range of MgCl₂ concentrations. 10X DreamTaq buffer includes 20 mM MgCl₂.

CERTIFICATE OF ANALYSIS

Tested for the absence of endo-, exodeoxyribonucleases, ribonucleases and functionally tested in amplification of 956 bp DNA fragment of a single-copy gene of human genomic DNA.

Quality authorized by:

 Jurgita Zilinskiene

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Please refer to www.thermoscientific.com/onebio for Material Safety Data Sheet of the product.

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Rev.7