

MegaBlock® 96 Well



Frei von humaner DNA,
DNase/RNase und Endotoxinen

MegaBlock® 96 Well

Der MegaBlock® 96 Well erfüllt alle wichtigen Anforderungen, wenn Proben bis zu 2,2 ml in automatisierten Systemen verarbeitet werden oder wenn große Mengen an Rückstellproben anfallen.

- Alphanummerische Beschriftung der Wells
- Hohe Sicherheit durch eine 100%ige Dichtheitsprüfung jedes Wells
- Ideal zur Langzeitlagerung von Proben
- Auch für Heatsealing-Systeme geeignet

MegaBlock® 0,5/1,2/2,2 ml, PP

- Lösungsmittelbeständig, auch gegen DMSO
- Autoklavierbar
- Quadratische Kavitäten mit Rundbodenwells
- Zur Rückstellung von pharmazeutischen Proben
- Für die DNA Isolierung, Enzymassays und Zellkulturanwendungen

MegaBlock® 1,2 ml, PS klar

- Leichte optische Kontrolle durch transparente Wells
- Optimal zur Langzeitlagerung von Blutproben
- Verschlussmatte durchstechbar

MegaBlock®

	Bestell-Nummer	Ausführung	Volumen/Material	Verpackung
NEU	82.1969.002	runde, erhabene Wells	0,5 ml PP	56 Stück / Karton
	82.1970.002	runde, erhabene Wells	1,2 ml PS klar	32 Stück / Karton
	82.1971.002	runde, erhabene Wells	1,2 ml PP	32 Stück / Karton
	82.1972.002	quadratische Wells	2,2 ml PP	32 Stück / Karton

Zur Abdeckung stehen Folien und Matten zur Verfügung

SARSTEDT AG & Co.
Postfach 12 20 · D-51582 Nümbrecht
Telefon (+49) 0 22 93 30 50
Telefax (+49) 0 22 93 305-282
☎ Service 0800 (Deutschland)
Telefon (0800) 0 83 30 50
info@sarstedt.com
www.sarstedt.com

 **SARSTEDT**

MegaBlock® 96 Well



Free from human DNA,
DNase/RNase, endotoxins

MegaBlock® 96 Well

The MegaBlock® 96 Well is ideal for automated processing of sample volumes up to 2.2 ml or handling large quantities of backup samples.

- Wells with alphanumeric coordination system
- Each well 100% leak-proof tested
- Ideal for long-term storage of samples
- Also suited for heat sealing

MegaBlock® 0.5/1.2/2.2 ml, PP

- Solvent-resistant, even to DMSO
- Autoclavable
- Square cavities with round bottom wells
- For backup storage of pharmaceutical substances
- For DNA isolation, enzyme and tissue culture applications

MegaBlock® 1.2 ml, PS transparent

- Transparent for viewing
- Ideal for long-term storage of blood samples
- Pierceable mat cover

MegaBlock®

	Order No.	Description	Volume/Material	Packaging
NEW	82.1969.002	round, raised wells	0.5 ml PP	56 pcs./case
	82.1970.002	round, raised wells	1.2 ml PS transparent	32 pcs./case
	82.1971.002	round, raised wells	1.2 ml PP	32 pcs./case
	82.1972.002	square wells	2.2 ml PP	32 pcs./case

Sealing foils/mats available

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305-0
Fax +49 2293 305-122
info@sarstedt.com
www.sarstedt.com