

5 ml Micro Tubes

The intermediate volume option for labs

NEW!



The Micro Tube for medium volumes:

- Convenient handling of contents between 2 and 5 ml
- Reliable SafeSeal-closure
- PCR Performance quality for more challenging applications



 **SARSTEDT**

5 ml Micro Tubes

Reliable handling of filling volumes between 2 and 5 ml

Routine applications in the lab like dilution series, cell breakdown or sample backup frequently require filling volumes above 2 ml. Due to their slender design or screw caps, larger volume tubes are sometimes not ideally suited for these purposes.

The new Sarstedt Micro Tubes combine the advantages of the SafeSeal Micro Tubes and 5 ml screw cap tubes in a single product. Easy to open with just one hand thanks to the ergonomic design, the snap-on lids with lip seal and retaining cams provide safe sealing characteristics even under thermal strain. The wide hinge enables easy and precise closure. Mixing up separate screw caps is excluded.



Advantages of the 5 ml Micro Tube at a glance:

- Robust polypropylene – centrifugation to 21,000 x g, autoclavable
- Moulded graduation – easy volume control
- Writing space and frosted lid – safe marking and quick identification
- Low height – minimal contamination risk as the pipette cone does not immerse into the tube
- Compact size – space-saving backup storage of samples, media or cultures

Ordering Information

Order No.	Description	Packaging
72.701	5 ml Micro Tube	100 / bag 200 / case
72.701.400	5 ml Micro Tube, PCR Performance Tested	100 / bag 200 / case



Accessories

Order No.	Description	Format	External dimensions (mm)	Packaging
93.850.001	Polycarbonate rack for 18 pcs. 5 ml Micro Tubes	6 x 3	137 x 70 x 40 (L x W x H)	1 / case
93.1431.001	Polycarbonate rack for 48 pcs. 5 ml Micro Tubes	12 x 4	257 x 90 x 40 (L x W x H)	1 / case

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305-0
Fax +49 2293 305-3992
export@sarstedt.com
www.sarstedt.com

