

REGANESTH® Spinal Needles

with ISO 80369-6 NRFit™ Non-Luer Connectors

NRFit™

Non-Luer Connector in accordance with ISO 80369-6




Superior manufacturing and design combined with patient safety

- Ergonomically designed needle hub
- CSF chamber with magnifying effect
- ISO 80369-6: new requirements for small bore connectors in the field of neuraxial applications and peripheral nerve blocks
- Reducing the risks of accidental misconnections of different supply lines to different access routes

REGANESTH® Spinal Needles

The Needle Hub

» unchanged cannula design




NRFit™
Non-Luer Connector in accordance with ISO 80369-6

- » optimised hub
- » comfortable handling
- » improved grip
- » brilliant CSF chamber

Needle Precision made by SARSTEDT

Pencil-Point



Quincke



Pencil-Point Ordering Information

Order No.	Description	Colour	Packaging
79.8081.514	REGANESTH® Spinal needle pencil-point NRFit 27G x 90 mm +Int.	●	20
79.8081.524	REGANESTH® Spinal needle pencil-point NRFit 27G x 38 mm	●	20
79.8081.534	REGANESTH® Spinal needle pencil-point NRFit 27G x 103 mm +Int.	●	20
79.8081.544	REGANESTH® Spinal needle pencil-point NRFit 27G x 120 mm +Int.	●	20
79.8081.804	REGANESTH® Spinal needle pencil-point NRFit 25G x 90 mm +Int.	●	20
79.8081.434	REGANESTH® Spinal needle pencil-point NRFit 25G x 103 mm +Int.	●	20
79.8081.844	REGANESTH® Spinal needle pencil-point NRFit 25G x 120 mm +Int.	●	20
79.8081.814	REGANESTH® Spinal needle pencil-point NRFit 24G x 90 mm +Int.	●	20
79.8081.824	REGANESTH® Spinal needle pencil-point NRFit 24G x 120 mm +Int.	●	20
79.8081.834	REGANESTH® Spinal needle pencil-point NRFit 24G x 150 mm +Int.	●	20
79.8081.214	REGANESTH® Spinal needle pencil-point NRFit 22G x 90 mm	●	20
79.8081.224	REGANESTH® Spinal needle pencil-point NRFit 22G x 120 mm	●	20

Quincke Ordering Information

Order No.	Description	Colour	Packaging
79.8082.254	REGANESTH® Spinal needle Quincke NRFit 25G x 90 mm +Int.	●	20
79.8082.124	REGANESTH® Spinal needle Quincke NRFit 22G x 120 mm	●	20
79.8082.184	REGANESTH® Spinal needle Quincke NRFit 22G x 180 mm	●	20
79.8082.024	REGANESTH® Spinal needle Quincke NRFit 22G x 38 mm	●	20
79.8082.214	REGANESTH® Spinal needle Quincke NRFit 22G x 50 mm	●	20
79.8082.224	REGANESTH® Spinal needle Quincke NRFit 22G x 90 mm	●	20
79.8082.204	REGANESTH® Spinal needle Quincke NRFit 20G x 90 mm	●	20

NRFit™ is a trademark of GEDSA and is used with their permission.

Distributed by:

SARSTEDT AG & Co. KG
P.O. Box 12 20
D-51582 Nümbrecht
Phone: +49 2293 305 0
Fax: +49 2293 305 3992
export@sarstedt.com
www.sarstedt.com

