

Magnetic Stand

Catalog # (P/N): AM10026
One Unit Includes: One single place magnetic stand

USER INFORMATION

Product Description: The Ambion® Single Place Magnetic Stand accommodates one 1.5 mL microcentrifuge tube. With this stand, centrifugation is replaced by the attraction between magnetic capture beads in a solution and the magnet in the stand, resulting in the beads and their captured biomolecules collecting on one side of the microcentrifuge tube.

Applications: The Magnetic Stand is intended for use when processing low numbers of samples in magnetic bead-based kits, such as MICROBExpress™ Bacterial mRNA Purification Kit (P/N AM1905) the MICROBEnrich™ Kit (P/N AM1901), or MagMAX™ Viral Isolation Kits (P/N AM1929 and AM1939).

QUALITY CONTROL

Functional Testing: The Magnetic Stand is functionally tested for magnetic separation of beads in a 1.5 mL microcentrifuge tube.

OTHER INFORMATION

Material Safety Data Sheets: Material Safety Data Sheets (MSDSs) can be printed or downloaded from product-specific links on our website at the following address: www.ambion.com/techlib/msds. Alternatively, e-mail your request to MSDS_Inquiry_CCRM@appliedbiosystems.com. Specify the catalog or part number(s) of the product(s), and we will e-mail the associated MSDSs unless you specify a preference for fax delivery. For customers without access to the internet or fax, our technical service department can fulfill MSDS requests placed by telephone or postal mail. (Requests for postal delivery require 1–2 weeks for processing.)

Warranty and Liability: *For research use only. Not for use in diagnostic procedures.*

Applied Biosystems is committed to delivering superior product quality and performance, supported by industry-leading global service and technical support teams. Warranty information for the accompanying consumable product is available at www.ambion.com/info/warranty in "Limited Warranty for Consumables," which is subject to the exclusions, conditions, exceptions, and limitations set forth under the caption "EXCLUSIONS, CONDITIONS, EXCEPTIONS, AND LIMITATIONS" in the full warranty statement. Please contact Applied Biosystems if you have any questions about our warranties or would like information about post-warranty support.

Information in this document is subject to change without notice. Applied Biosystems assumes no responsibility for any errors that may appear in this document.

Applied Biosystems disclaims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall Applied Biosystems be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

Applied Biosystems, AB (Design), and Ambion are registered trademarks, and MagMAX, MICROBEnrich, and MICROBExpress, are trademarks of Applied Biosystems Corporation or its subsidiaries in the US and/or certain other countries. All other trademarks are the sole property of their respective owners.

Trademarks, Patents, and Licensing:

© 2008 Ambion, Inc. All rights reserved. 4386454A