

# Transfusion devices



- For the transfusion of blood products
- Use of fine-pored sieve filters for effective elimination of micro aggregates
- High filter capacity and high flow rate due to a large filter surface





- 1 TRANSMED 10** – Integrated cascading filter 175 µm, 40 µm and 10 µm
- 2 TRANSMED 40** – Integrated cascading filter 175 µm and 40 µm
- 3 TRANSFLOW** – Integrated fine filter with a pore size of 40 µm
- 4 TRANSFLOW-N** – Integrated fine filter with a pore size of 40 µm

## Ordering information

Order no.	Description	Pck.
79.8200.000	<b>1</b> TRANSMED 10	100
79.8200.100	<b>2</b> TRANSMED 40	100
79.8200.440	<b>3</b> TRANSFLOW	100
79.8200.490	<b>4</b> TRANSFLOW-N	100

### Distributed by:

SARSTEDT AG & Co. KG  
 P.O. Box 12 20  
 D-51582 Nümbrecht  
 Phone: +49 2293 305 0  
 Fax: +49 2293 305 3992  
 export@sarstedt.com  
 www.sarstedt.com