

## BRC-5300T High Speed Bench-Top Touch Panel Centrifuge Rotor Types



Order No.	Rotor Type	Max. Speed (rpm)	Max. RCF (xg)	Max. Capacity
INO53T-1	Swing Rotor (Round Bucket)	4000	3180	4x300ml
INO53T-2	Swing Rotor (Square Bucket)	4000	3150	4x300ml
INO53T-3	Swing Rotor	4000	3150	4x24x7ml
INO53T-4	Swing Rotor	4000	2750	4x24x5ml
INO53T-5	Swing Rotor	4000	3150	4x16x10ml
INO53T-6	Swing Rotor	4000	3410	4x22x7ml
INO53T-7	Swing Rotor	4000	2840	4x22x5ml
INO53T-8	Swing Rotor	4000	3140	4x18x10ml
INO53T-9	Microplate Rotor (Box Type)	4000	2770	2x4x96 well
INO53T-10	Microplate Rotor (Fame Type)	4000	2910	2x4x96 well
INO53T-11	Microplate Rotor (Fame Type)	4000	2300	2x3x48 well
INO53T-12	Swing Rotor	5000	4730	4x1x50ml
INO53T-13	Swing Rotor	5000	4730	4x1x100ml
INO53T-14	Swing Rotor	4000	3020	4x2x100ml
INO53T-15	Swing Rotor	4000	3200	4x2x50ml
INO53T-16	Swing Rotor	4000	3200	4x4x15ml
INO53T-17	Swing Rotor	4000	3200	4x6x15ml
INO53T-18	Swing Rotor	4000	3200	4x8x15ml
INO53T-19	Swing Rotor	4000	2770	4x8x10/7ml
INO53T-20	Swing Rotor	4000	2380	4x8x5ml adapter
INO53T-21	Swing Rotor	4000	2770	4x6x10/7ml
INO53T-22	Swing Rotor	4000	2380	4x6x5ml adapter
INO53T-23	Swing Rotor	4000	2770	4x4x10/7ml
INO53T-24	Swing Rotor	4000	2380	4x4x5ml adapter
INO53T-25	Swing Rotor	4000	2520	4x16x5ml
INO53T-26	Swing Rotor	4000	2880	4x16x7ml
INO53T-27	Swing Rotor	4000	2880	4x12x10ml
INO53T-28	Fixed Rotor	14000	12070	4x8 PCR
INO53T-29	Fixed Rotor	13000	16080	6x8 PCR
INO53T-30	Fixed Rotor	14000	13390	8x8 PCR
INO53T-31	Fixed Rotor	13000	17220	12x8 PCR
INO53T-32	Fixed Rotor	14000	19960	40x0.5ml
INO53T-33	Fixed Rotor	18500	23990	12x1.5ml
INO53T-34	Fixed Rotor	15000	20600	24x1.5ml
INO53T-35	Fixed Rotor	13500	19340	30x1.5ml
INO53T-36	Fixed Rotor	11000	12840	48x1.5ml
INO53T-37	Fixed Rotor	15000	19350	16x5ml
INO53T-38	Fixed Rotor	14000	16460	6x10ml
INO53T-39	Fixed Rotor	13000	17030	12x10ml
INO53T-40	Fixed Rotor	11000	13950	16x10ml
INO53T-41	Fixed Rotor	11000	12730	8x15ml
INO53T-42	Fixed Rotor	9000	9590	12x15ml
INO53T-43	Fixed Rotor	5000	3500	24x15ml
INO53T-44	Fixed Rotor	5000	3830	30x15ml
INO53T-45	Fixed Rotor	12000	14000	6x30ml
INO53T-46	Fixed Rotor	12000	16050	6x50ml (Conical)
INO53T-47	Fixed Rotor	12000	15960	6x50ml (Round)
INO53T-48	Fixed Rotor	5000	3860	12x50ml
INO53T-49	Fixed Rotor	5000	2910	24x50ml
INO53T-50	Fixed Rotor	11000	13560	4x85ml
INO53T-51	Fixed Rotor	10000	10310	4x100ml
INO53T-52	Fixed Rotor	8000	7280	6x100ml
INO53T-53	Fixed Rotor	4000	2970	12x100ml
INO53T-54	Swing Rotor	12000	15800	24 Capillaries
INO53T-55	Vertical Rotor	15000	19920	4x5ml
INO53T-56	Vertical Rotor	16000	16450	16x5ml
INO53T-57	Vertical Rotor	14000	19750	8x30ml