

BRC-5300UT High Speed Bench-Top Touch Panel Centrifuge Rotor Types



Order No.	Rotor Type	Max. Speed (rpm)	Max. RCF (xg)	Max. Capacity
INO53UT-1	Swing Rotor (Round Bucket)	4000	3180	4x300ml
INO53UT-2	Swing Rotor (Square Bucket)	4000	3150	4x300ml
INO53UT-3	Swing Rotor	4000	3150	4x24x7ml
INO53UT-4	Swing Rotor	4000	2750	4x24x5ml
INO53UT-5	Swing Rotor	4000	3150	4x16x10ml
INO53UT-6	Microplate Rotor (Fame Type)	4000	2910	2x4x96 well
INO53UT-7	Microplate Rotor (Box Type)	4000	2770	2x4x96 well
INO53UT-8	Microplate Rotor	4000	2300	2x3x48 well
INO53UT-9	Swing Rotor	4000	3410	4x22x7ml
INO53UT-10	Swing Rotor	4000	2840	4x22x5ml
INO53UT-11	Swing Rotor	4000	3140	4x18x10ml
INO53UT-12	Swing Rotor	4000	4730	4x1x50ml
INO53UT-13	Swing Rotor	5000	4730	4x1x100ml
INO53UT-14	Swing Rotor	5000	3020	4x2x100ml
INO53UT-15	Swing Rotor	4000	3200	4x2x50ml
INO53UT-16	Swing Rotor	4000	3200	4x4x15ml
INO53UT-17	Swing Rotor	4000	3200	4x6x15ml
INO53UT-18	Swing Rotor	4000	3200	4x8x15ml
INO53UT-19	Swing Rotor	4000	2770	4x8x10/7ml Vacuum Tube
INO53UT-20	Swing Rotor	4000	2380	4x8x5ml Vacuum Tube
INO53UT-21	Swing Rotor	4000	2770	4x6x10/7ml Vacuum Tube
INO53UT-22	Swing Rotor	4000	2380	4x6x5ml Vacuum Tube
INO53UT-23	Swing Rotor	4000	2770	4x4x10/7ml Vacuum Tube
INO53UT-24	Swing Rotor	4000	2380	4x4x5ml Vacuum Tube
INO53UT-25	Swing Rotor	4000	2520	4x16x5ml Vacuum Tube
INO53UT-26	Swing Rotor	4000	2880	4x16x7ml Vacuum Tube
INO53UT-27	Swing Rotor	4000	2880	4x16x10ml Vacuum Tube
INO53UT-28	Swing Rotor	4000	2880	4x16x10ml Vacuum Tube
INO53UT-29	Fixed Rotor	16000	15760	4x8 PCR
INO53UT-30	Fixed Rotor	15000	21420	6x8 PCR
INO53UT-31	Fixed Rotor	16000	17480	8x8 PCR
INO53UT-32	Fixed Rotor	15000	22930	12x8 PCR
INO53UT-33	Fixed Rotor	15000	22920	40x0.5ml
INO53UT-34	Fixed Rotor	21000	30910	12x1.5ml
INO53UT-35	Fixed Rotor	16000	23440	24x1.5ml
INO53UT-36	Fixed Rotor	14000	20800	30x1.5ml
INO53UT-37	Fixed Rotor	13000	17930	48x1.5ml
INO53UT-38	Fixed Rotor	16000	22020	16x5ml
INO53UT-39	Fixed Rotor	16000	21500	6x10ml
INO53UT-40	Fixed Rotor	15000	22680	12x10ml
INO53UT-41	Fixed Rotor	13000	19490	16x10ml
INO53UT-42	Fixed Rotor	13000	17790	8x15ml
INO53UT-43	Fixed Rotor	11000	14330	12x15ml
INO53UT-44	Fixed Rotor	5000	3500	24x15ml
INO53UT-45	Fixed Rotor	5000	3830	30x15ml
INO53UT-46	Fixed Rotor	14000	19060	6x30ml
INO53UT-47	Fixed Rotor	13000	18840	6x50ml (Conical)
INO53UT-48	Fixed Rotor	13000	18730	6x50ml (Round)
INO53UT-49	Fixed Rotor	5000	3860	12x50ml
INO53UT-50	Fixed Rotor	4000	2970	24x50ml
INO53UT-51	Fixed Rotor	13000	18940	4x85ml
INO53UT-52	Fixed Rotor	12000	14850	4x100ml
INO53UT-53	Fixed Rotor	10000	11380	6x100ml
INO53UT-54	Fixed Rotor	4000	2970	12x100ml
INO53UT-55	Fixed Rotor	12000	15800	24 Capillaries
INO53UT-56	Swing Rotor	15000	19920	4x5ml
INO53UT-57	Vertical Rotor	16000	16450	16x5ml
INO53UT-58	Vertical Rotor	14000	19750	8x30ml