

Universal Container



Die bruch sicheren, besonders dichten 25 ml-Röhren sind ideal geeignet für jeden Verarbeitungsschritt von analytischem und diagnostischem Material – von der Abfüllung über den Transport, die Vorbereitung und Analyse bis hin zur Lagerung.

Hergestellt aus hochtransparentem, robustem Polypropylen und mit praktischem Stehfuß garantieren sie bei jedem Arbeitsschritt eine sichere Inspektion und Handhabung der Probe.

Produktvorteile auf einen Blick

- Transparentes Polypropylen – zuverlässige visuelle Kontrolle
- Robustes Material – zentrifugierbar bis 4.500 x g
- Optimale Transportsicherheit – Dichtigkeit entspricht der Verpackungsvorschrift P 650 der ADR und IATA
- Montierter Schraubverschluss – praktische, hygienische Entnahme
- Optionales Papieretikett mit Schriftfeld und Graduierung – sichere Beschriftung und schnelle Volumenabschätzung

Bestellinformation 25ml Röhre 90 x 25mm

Bestell-Nr.	Verschlussfarbe	Etikett	steril	Material	Verpackung
60.9922.241				PP	500/Beutel
63.9922.252		✓		PP	500/Beutel
60.9922.243			✓	PP	500/Beutel
63.9922.254		✓	✓	PP	500/Beutel
63.9922.248		✓	✓	PP	1/Blister 300/Beutel

SARSTEDT AG & Co. KG
 Postfach 12 20 · D-51582 Nümbrecht
 Telefon +49 22 93 305 - 0
 Telefax +49 22 93 305 - 3450
Service 0800 (Deutschland)
 Telefon (0800) 0 83 305 - 0
 info@sarstedt.com
 www.sarstedt.com

Universal Container



The break-resistant, especially leak-proof 25 ml tubes are ideally suited for the complete processing workflow of analytical and diagnostic material – from sample filling to transport, preparation and analysis to storage.

Produced from highly transparent, durable polypropylene and thanks to a convenient skirted base, safe inspection and handling of the sample is guaranteed at every work step.

Product benefits at a glance

- Transparent polypropylene – reliable visual control
- Durable material – centrifugation up to 4,500 x g
- Optimal transport safety – leakproof in accordance with P 650 packaging instructions of the ADR and IATA
- Assembled screw cap – convenient, hygienic removal
- Optional paper label with writing space and graduation – safe marking and quick volume estimation

Ordering Information 25 ml Tube 90 x 25mm

Order No.	Cap colour	Label	Sterile	Material	Packaging
60.9922.241				PP	500/bag
63.9922.252		✓		PP	500/bag
60.9922.243			✓	PP	500/bag
63.9922.254		✓	✓	PP	500/bag
63.9922.248		✓	✓	PP	1/blister 300/bag

SARSTEDT AG & Co. KG
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305-0
Fax +49 2293 305-3992
export@sarstedt.com
www.sarstedt.com