

# High-quality water by design

Take advantage of special pricing for 20 L and 200 L Gibco WFI for Cell Culture

## Gibco™ Water for Injection (WFI) for Cell Culture

For further manufacturing

- Multi-compendial
  - Meets standards for packaged sterile purified water and sterile water for injection for further manufacturing
  - Compliant with United States Pharmacopeia (USP) and European Pharmacopoeia (EP) specifications
- Manufactured in ISO-certified, cGMP, and FDA-registered facilities
- Available: cGMP warehousing and just-in-time delivery services


### Promotional pricing only for these Gibco WFI for Cell Culture sizes

Size	Promotional price	List price	Cat. No.
20 L	\$250	\$408	A1287305
200 L	\$750	\$1,638	A1287306



## 20 L packaging design

(Thermo Scientific™ CX5-14 film)




### Materials used in 20 L standard Thermo Scientific™ BPC

Line	Tubing (ID x OD x wall, total length)	Connector(s)
Line 1	Thermo Fisher Scientific fill line	Sealed after filling
Line 2	C-Flex 374 tubing 3/8 x 5/8 x 1/8 in., total length 5 feet (9.5 x 15.9 x 3.2 mm, total length 1.2 meters)	Polycarbonate 3/8 in. (9.5 mm) MPC series quick-connect body (female)
Line 3	C-Flex 374 tubing 1/4 x 3/8 x 1/16 in., total length 5 feet (6.4 x 9.5 x 1.6 mm, total length 1.2 meters)	Polycarbonate 1/4 in. (6.4 mm) MPC series quick-connect insert (male)

## 200 L packaging design

(Thermo Scientific CX5-14 film)



### Materials used in standard top-drain 200 L Thermo Scientific BPC

Line	Tubing (ID x OD x wall, total length)	Connector(s)
Line 1	C-Flex 374 tubing 3/8 x 5/8 x 1/8 in., total length 2 feet (9.5 x 15.9 x 3.2 mm, total length 0.61 meters)	Polycarbonate 3/8 in. (9.5 mm) MPC series quick-connect body (female)
Line 2	C-Flex 374 tubing 1/2 x 3/4 x 1/8 in., total length 2.5 feet (12.7 x 19.1 X 3.2 mm, total length 0.76 meters)	Polycarbonate 1/2 in. (12.7 mm) MPX series quick-connect insert (male)
Line 3	Thermo Fisher Scientific fill line	Sealed after filling

Contact your sales representative to take advantage of this promotion.

Find out more at [thermofisher.com/wfi](https://thermofisher.com/wfi)

**ThermoFisher**  
SCIENTIFIC

**For Research Use or Further Manufacturing. Not for diagnostic use or direct administration into humans or animals.**

© 2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. C-Flex is a trademark of Saint-Gobain Performance Plastics Corporation. **COL18706 0819**

\* This promotion is open to customers in North, Central, and South America and Europe. Qualifying orders must be placed through your sales representative. Discount will apply to qualifying orders received by Life Technologies no later than December 31, 2020, or until promotional supplies are depleted, whichever comes first. Cannot be combined with other discounts or promotions. Offer void where prohibited, licensed, or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply.