

Labor-Technik-Göttinger

## **Blockthermostat**

**Quality – Made in Germany** 



## Technical data:

- · Dimensions thermostat: 290 mm x 175 mm x 110 mm  $(W \times D \times H)$
- · Power supply: 100 - 240 V AC, 50/60 Hz, 60 W
- Temperature control between 27 °C and 45 °C
- Accuracy: ± 0.2 °C
- Temperature distribution: ± 0.2 °C
- Heating up time from 22 °C to 37 °C: approx. 40 min.
- · Safety class II

## **Designed to maintain the** temperature of biological materials

## Features:

- Temperature adjustable between 27 °C and 45 °C
- Homogeneous temperature distribution
- Heating up time only approx. 40 min.
- Use of different heating blocks

**REF** 10719 (without heating block)

13855 (with heating block: REF REF 10264)

REF 13896 (with heating block: REF 10277)

CE