

V-Monovette® Urine

The vacuum system for enclosed urine transfer

NEW!



V-Monovette® Urine Product Range

Comprehensive product range

The V-Monovette® Urine is suited for a variety of applications and available in the size options 100 x 15, 75 x 13 mm and 100 x 13 mm with rounded base and 100 x 15 mm with conical base.

For urine collection, the system includes a 100 ml urine cup and a 3 l urine container with integrated sampling port for enclosed urine transfer to the V-Monovette®. A transfer device for urine sampling from open containers is also available.

Optimal hygiene and convenient handling

An enclosed system significantly improves hygiene and is user friendly. With this system, it is no longer necessary to open a urine cup or container to fill the tube so that contamination with urine is practically excluded.



Ordering Information

	Order No.	Description	Length/ø in mm	Packaging pcs. inner case/box
Urinalysis	11.2252.001	V-Monovette® Urine 4 ml*	75 x 13	50/500
	11.2352.001	V-Monovette® Urine 6 ml*, rounded base	100 x 13	50/500
	11.2452.001	V-Monovette® Urine 10 ml*, rounded base	100 x 15	50/500
	11.2552.001	V-Monovette® Urine 9.5 ml*, conical base	100 x 15	50/500
Urine Culture	11.2253.001	V-Monovette® Urine 4 ml with stabilizer*, rounded base	75 x 13	50/500
	11.2453.001	V-Monovette® Urine 10 ml with stabilizer*, rounded base	100 x 15	50/500
Accessory	11.1240	Transfer device		50/500
	75.562.400	100 ml Cup with integrated transfer device, assembled lid, sterile, with warranted sterility due to safety label		200/box
	77.575.400	3 l Urine container with integrated transfer device		30/box
Accessory Urine Collection from Catheter	14.1208.700	Needle holder		100/1.000
	85.1638.935	Safety-Multifly® Needle with assembled tubeholder, 21 G x 1¼" tubing 8"		100/500

* All versions of the V-Monovette® Urine are sterile

SARSTEDT Inc.
1025 St. James Church Road
P.O. Box 468
Newton, NC 28658-0468
Tel. (800) 257-5101
Fax (828) 465-4003
sarstedt@bellsouth.net
www.sarstedt.com

