


S-Monovette® Rack



S-Monovette® Rack - A must for any laboratory!

- Available in three different sizes:
 - Rack D17, 50 apertures, Ø 17 mm
 - Rack D12, 50 apertures, Ø 12 mm
 - Rack D13, 50 apertures, Ø 13 mm
 - D13 and D17 also available **with handles**
- Especially suitable for space-saving sample backup storage
- Can be compiled to form sample distribution blocks
- Highly impact-resistant, stackable, autoclavable
- Designed to accommodate all S-Monovette® options
- Click on our “Tube Finder“ at www.sarstedt.com for further information on compatible racks for tubes of various sizes.

Available in 9 colours!

Rack D12 (L x W x H) 152 x 80 x 45 mm
50 apertures, e.g. suitable for S-Monovettes and tubes from
Ø 8 mm to Ø 12 mm

Colour	Order No.	Packaging
grey	93.854	20/box
red	93.854.121	20/box
yellow	93.854.122	20/box
green	93.854.123	20/box
blue	93.854.124	20/box
orange	93.854.125	20/box
grey	93.854.200 double block rack = 100 apertures connected at lateral ends*	12/box



Rack D13 (L x W x H) 180 x 96 x 45 mm
50 apertures, e.g. suitable for S-Monovettes and tubes from
Ø 8 mm to Ø 13 mm


Colour	Handles	Order No.	Packaging
grau	without	93.853	20/box
	with	93.855	24/box
rot	without	93.853.131	20/box
	with	93.855.131	24/box
gelb	without	93.853.132	20/box
	with	93.855.132	24/box
grün	without	93.853.133	20/box
	with	93.855.133	24/box
blau	without	93.853.134	20/box
	with	93.855.134	24/box
orange	without	93.853.135	20/box
	with	93.855.135	24/box
grey (other colours available on request)	without	93.853.200 double block rack = 100 apertures connected at lateral ends*	10/box
grey (other colours available on request)	without	93.853.300 double block rack = 100 apertures connected at front ends*	10/box

* Other colours available on request.

At a minimum quantity of 100 pcs., the ends can be printed to your specifications (except for options with handles).

Double block rack

Two D17 or D13 racks are firmly connected either at their front or lateral ends to form a rack with 100 apertures. Particularly convenient for sample backup storage in refrigerators.



Rack D17 (L x W x H) 209 x 109 x 45 mm
50 apertures, e.g. suitable for S-Monovettes and tubes from
Ø 13 mm to Ø 17 mm

Colour	Handles	Order No.	Packaging
grey	without	93.852	20/box
	with	93.851	28/box
red	without	93.852.171	20/box
	with	93.851.171	28/box
yellow	without	93.852.172	20/box
	with	93.851.172	28/box
green	without	93.852.173	20/box
	with	93.851.173	28/box
blue	without	93.852.174	20/box
	with	93.851.174	28/box
orange	without	93.852.175	20/box
	with	93.851.175	28/box
brown	without	93.852.168	20/box
violet	without	93.852.169	20/box
white	without	93.852.170	20/box
grey (other colours available on request)	without	93.852.200 double block rack = 100 apertures connected at lateral ends*	16/box
grey (other colours available on request)	without	93.852.300 double block rack = 100 apertures connected at front ends*	12/box
grey (other colours available on request)	with	93.851.200 double block rack = 100 apertures connected at lateral ends*	14/box

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 3050
Fax +49 2293 305-122
info@sarstedt.com
www.sarstedt.com