

vivantis

DNA AMPLIFICATION PRODUCT Product Datasheet

**10X
ViBuffer S**

Product No : RB0203
Packaging : 5 x 1ml

Lot :
Expiry Date :

Store at -20°C info@vivantechnologies.com

Description:

General reagents used for various PCR reactions.

Application:

Generally used in various PCR reactions with Vivantis thermostable DNA Polymerases.

For Research Use Only



Buffer Composition:

160mM (NH₄)₂SO₄,
500mM Tris-HCl (pH 9.2 at 22°C),
17.5mM MgCl₂ and 0.1% Triton X-100.

The buffer is optimized for use with 0.35mM of each dNTP

Quality Assurance:

Functionally tested in PCR with *Taq* DNA Polymerase and *Pfu* DNA Polymerase.

vivantis | www.vivantechnologies.com

DSRB0203_rev0_010311