

Aspirator 3

Quality – Made in Germany



Precise aspiration pump for follicular puncture

Basic equipment:

- Aspirator 3
- Foot pedal
- · Mains adapter
- 2 x overflow bottle, 2 x blocking filter and tubing

REF 14193 (without hotplate)

REF 13907 (with integrated hotplate)

Advantages of Aspirator 3:

- Approved medical device
- Easy handling by electrical pedal
- Long durability
- High reliability
- Very silent
- Constant control of negative pressure
- Safe by overflow protection

Picture shows device with integrated hotplate (optional), heating block and puncture set from Labotect

Optional accessory:

Integrated hotplate

REF | 13907 (Aspirator 3 with hotplate)

Movable support stand



Technical data:

- According to Directive 93/42/EEC: Class IIa
- EN 10079 conform: high vacuum, low flow
- · Volume of overflow vessel: 100 ml
- Dimensions:
- $\sim 408 \times 130 \times 250 \text{ (w x h x d in mm)}$
- Weight: ~ 7.6 kg
- Connection for aspiration tubing: Ø 2 to 4 mm
- Power supply: 230 V AC, 50/60 Hz or 115 V AC, 50/60 Hz, 85 VA
- Safety class: I
- Device type: BF
- Hotplate (optional):
 - Temperature range: 30 °C 45 °C
 - Accuracy: ± 0.2 °C
 - Temperature distribution: ± 0.5 °C
 - Heating up time: < 10 min.
- Connection for potential compensation



(€0123

Rev. 5_07/2015