



**Benchtap
Centrifuge**

EBA 20

Practical and handy, the EBA 20 is the ideal centrifuge for small sample volume.

Complete with an 8-place angle rotor.

Microprocessor-controlled.

Max. speed: 6,000 min⁻¹

Max. RCF: 3,421

Variable – it accommodates various tubes as well as glass tubes without adapters.



Ergonomically arranged controls and displays.

Entry of the parameters via foil keypad.

Speed RPM: Input in increments of 100.

Time: Input in minutes, max. 99 min.

Impulse key: For short centrifuging operations.

Display of the actual values during centrifugation.

“Open lid“ key



8 x 15 ml standard tubes
(17 x 100 mm)



8 x 8-9 ml Greiner Vacuettes or
8 x 8-10 ml Vacutainers (16 x 100 mm)



4 x 15 ml Falcon tubes
(17 x 120 mm)



4 x 8 ml Vacutainers
(16 x 125 mm)



8 x 9-10 ml Sarstedt
Monovettes (16 x 92 mm)

Adapters for a variety of additional blood collection tubes on request.

Safety

Lid locking

Lid holding

Emergency lid lock release

Motor overheating protection

Transparent lid

Conformity with safety regulations like IEC 1010-2-2, the relevant EC standards, and the DIN norms.

Technical Details

Power supply*:	208-240 V 1~
Frequency:	50-60 Hz
Consumption:	65 VA
Radio interference:	EN 55011 ISM class B
Immunity:	as per EN 50082 - 1
Radius/max. speed/max. RCF:	85 mm/6,000 min ⁻¹ /3,421
Running time:	1-99 min., cont. run, short cycle mode
Dimensions in mm (HxWxD):	216 x 231 x 292
Weight:	approx. 4 kg

*other voltages on request

Ordering Information

EBA 20

Benchtop centrifuge with rotor
Cat.-No. 2002

Design:

Housing and lid are made of scratch and impact resistant plastic. Angle rotor of steel with buckets of glass fibre reinforced polyamide.

Hettich
ZENTRIFUGEN

Andreas Hettich GmbH & Co. KG
Gartenstr. 100 · D-78532 Tuttlingen / Germany
Phone +49 (0) 7461 705-0 · Fax +49 (0) 7461 705-122

International Sales +49 (0) 7461 705-201

International Service +49 (0) 7461 705-203

www.hettichlab.com · info@hettichlab.com · service@hettichlab.com

2002.0404 © by Hettich Tuttlingen · Printed in Germany · Design, form and colour modifications reserved