

# Hot Plate 062

**Quality – Made in Germany** 



## Hot Plate for laboratory and scientific use

#### Features:

- Temperature control between 27 °C and 45°C
- Homogeneous temperature distribution
- Planar surface for optimal heat transfer to heating blocks, culture dishes and slides
- · Short heating up time
- · Very compact und extremely flat
- Space saving due to separate control unit

#### Technical data:

- Dimensions: 274 x 214 x 18 mm (w x d x h)
- Weight: 2.4 kg
- Accuracy: ± 0.2 °C
- Temperature distribution: ± 0.3 °C
- Heating up time: < 10 min.</li>
- Power supply: 100 240 V AC, 50/60 Hz, 60 W
- Temperature adjustable between 27 °C and 45 °C
- · Safety class II

### Scope of supply:

- Hot plate
- · Controller with mains cable

REF | 13854

## Additional equipment:

Retractor frame

REF | 15943



Rev. 4\_01/2017