TaqMan™ RNase P Instrument Verification Plate, 96-well

Catalog Number 4310982

Pub. No. 4314333 Rev. H



WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from **thermofisher.com/support**.

Contents and storage

Contents	Amount	Storage			
TaqMan™ RNase P Instrument Verification Plate, 96-well shipped with	1	−25°C to −15°C			
Optical Compression Pad		Do not store in a frost-free freezer.			

RNase P plate layout

	1	2	3	4	5	6	7	8	9	10	11	12	
Α													
В													
С													
D	NTC	NTC	NTC	NTC			STND 1.25K			STND 2.5K		STND 2.5K	
Ε	STND 5K	STND 5K	STND 5K	STND 5K	STND 10K	STND 10K	STND 10K	STND 10K	STND 20K	STND 20K	STND 20K	STND 20K	
F													
G													
Н													
	Unknown 5K							Unknown 10K					

Related documentation

For detailed information on the plate set-up and analysis procedures, refer to the appropriate Sequence Detection System user guide.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



The information in this guide is subject to change without notice.

DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a trademark of Roche Molecular Systems, Inc., used under permission and license.

Thermo Fisher SCIENTIFIC