

Prep-n-Go™ Buffer

Publication Part Number 4466622 Rev. B Revision Date June 2012

Part Number	Part	Storage Conditions
4467079	Prep-n-Go™ Buffer (untreated paper substrate), 5 tubes, 0.5 mL each	Room temperature, protect from light
4471406	Prep-n-Go™ Buffer (buccal swab), 1 bottle, 80 mL	

Note: For safety and biohazard guidelines, refer to the “Safety” section in the user guide for your AmpF ϕ STR® Direct PCR Amplification Kit. For every chemical, read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product description

Prep-n-Go™ Buffer is optimized for use with AmpF ϕ STR® Direct PCR Amplification Kits (for example, the AmpF ϕ STR® Identifier® Direct PCR Amplification Kit, 200 reaction kit Part no. 4467831; 1000 reaction kit Part no. 4408580; and other and other Applied Biosystems® STR direct chemistries).

Prep-n-Go™ Buffer is optimized for the following sample types:

- Part no. 4467079 optimized for untreated paper substrate
- Part no. 4471406 optimized for buccal swab substrate

For instructions on using the Prep-n-Go™ Buffer, refer to the user guide for your AmpF ϕ STR® Direct PCR Amplification Kit.

Performance characteristics

The Prep-n-Go™ Buffer is developed and manufactured by Life Technologies Corporation. We test each lot of kits to ensure that the buffer performs according to specifications and is free from human DNA contamination.

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

© 2012 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners.

For support visit www.lifetechnologies.com/support

www.lifetechnologies.com

