



Comments for pBAD/His A
4102 nucleotides

- araBAD* promoter region: bases 4-276
- Initiation ATG: bases 319-321
- Polyhistidine tag: bases 331-348
- Xpress™ epitope: bases 388-411
- Enterokinase recognition site: bases 397-411
- Multiple cloning site: bases 430-470
- rrnB* transcription termination region: bases 553-710
- Ampicillin ORF: bases 989-1849
- pBR322 origin: bases 1994-2667
- AraC ORF: bases 4076-3198

**Sac I* and *Pvu II* are not present in version C.