

BIOLGIX[®]



**NEW MEMBERS OF
DEEP WELL PLATES!**

For COVID-19 Sample Testing


•SBS Footprint •Autoclavable •Automated DNA Extraction Equipment

96 Square Wells, V Bottom



**2.2ml & 0.5ml volumes,
2.2ml plates can be used with
96 tip combo**

96 Square Wells,
U Bottom, Profile Concave



**2.2ml volumes,
can be used with magnet set**

96 Square Wells, U Bottom



1.6ml & 2.2ml volumes

96 Round Wells, U Bottom



1.0ml volume

24 Square Wells, V Bottom




10ml volume

96 Deep Well Silicone
Sealing Mats



**Round & square wells,
punctured & non-punctured**

Order Information



CAT #	DESCRIPTION
02-8005	0.5ml deep well plates, square 96 wells, V bottom, match with KingFisher Flex system, without cap, non-sterile, 10 Pieces/Pack, 100 Pieces/Case
02-2010	1.0ml deep well plates, round 96 wells, U bottom, without cap, non-sterile, 24 Pieces/Pack, 96 Pieces/Case
02-1016	1.6ml deep well plates, square 96 wells, U bottom, without cap, non-sterile, 24 Pieces/Pack, 96 Pieces/Case
02-1022	2.2ml deep well plates, square 96 wells, U bottom, without cap, non-sterile, 24 Pieces/Pack, 96 Pieces/Case
02-6022	2.2ml deep well plates, square 96 wells, U bottom, profile concave, can be used with 8 magnet set, without cap, non-sterile, 24 Pieces/Pack, 96 Pieces/Case
02-7022	2.2ml deep well plates, square 96 wells, V bottom, match with KingFisher Flex system, can be used with 96 tip comb, without cap, non-sterile, 10 Pieces/Pack, 50 Pieces/Case
02-9010	10ml deep well plates, square 24 wells, V bottom, match with KingFisher Flex system, without cap, non-sterile, 10 Pieces/Pack, 50 Pieces/Case
02-0096	96 tip comb, non-Sterile, 2 Pieces/Pack, 50 Pieces/Case
02-0196	96 tip comb, sterile, 2 Pieces/Pack, 50 Pieces/Case
02-0008	Magnet set, 2 Pieces/Pack, 100 Pieces/Case
02-1096	96 square deep well silicone sealing mats, for 96 square deep well plates, 10 Pieces/Pack, 50 Pieces/Case
02-1196	96 square deep well silicone sealing mats, can be punctured, for 96 square deep well plates, 10 Pieces/Pack, 50 Pieces/Case
02-2196	96 round well silicone sealing mats, can be punctured, for 1.0ml V-bottom & U-bottom Plates, 10 Pieces/Pack, 50 Pieces/Case
02-2296	96 round well silicone sealing mats, can be punctured, for 2.0ml deep well plates, 10 Pieces/Pack, 50 Pieces/Case
02-3024	24 square deep well silicone sealing mats, for 24 square deep well plate 10ml, 10 Pieces/Pack, 50 Pieces/Case

*Sterilized, individual wrapped, and with lid type optional.

Biologix Group Limited

+86-531-67802668 Inquiry@BiologixGroup.com www.BiologixGroup.com