



Comments for pRSET/CFP, EmGFP and BFP
3600 nucleotides

T7 promoter/priming site: bases 9-28

6xHis tag: bases 101-118

T7 gene 10 leader: bases 122-154

Xpress™ epitope: bases 158-181

EK cleavage site: bases 167-181

Fluorescent Protein (CFP, EmGFP, BFP): bases 209-928

T7 reverse priming site: bases 987-1006

T7 transcription terminator: bases 948-1084

f1 origin: bases 1148-1603

bla promoter: bases 1635-1739

Ampicillin (bla) resistance gene: bases 1734-2594

pUC origin: bases 2739-3412