


PCR PRODUCTS-PLATES



Half-skirted plates are easy to insert and remove


Features:

- Made of medical grade polypropylene.
- 0.2ml, 96 wells.
- Ultra thin walls with uniform thickness for accelerated heat transfer, shorter cycle times, and more consistent results.
- Alphanumeric well labeling for convenient identification.
- DNase, RNase, and endotoxin free.
- Fit most standard thermocyclers.
- Can withstand centrifugation of 6,000 RCF.

Thermocyclers



Sealing Films



Biologix Group Limited

+86-531-67802668 Inquiry@BiologixGroup.com www.BiologixGroup.com