



**Comments for pDEST™24
6961 nucleotides**

T7 promoter: bases 25-44

attR1: bases 71-195

Chloramphenicol resistance gene (Cm^R): bases 304-963

ccdB gene: bases 1305-1610

attR2: bases 1651-1775

GST tag: bases 1783-2451

T7 transcription termination region: bases 2466-2594

bla promoter: bases 3082-3180

Ampicillin (*bla*) resistance gene: bases 3181-4041

pBR322 origin: bases 4186-4859

ROP ORF: bases 5230-5421 (C)

C=complementary strand