



**Comments for pTracer™-CMV2  
6210 nucleotides**

CMV promoter: bases 209-863

CMV forward priming site: bases 769-789

Multiple cloning site: bases 901-1010

BGH reverse priming site: bases 1022-1039

BGH polyadenylation signal: bases 1025-1223

EF-1a promoter: bases 1304-2487

EM-7 promoter: bases 2603-2658

GFP-Zeocin™ resistance gene: bases 2677-3753

GFP ORF: bases 2677-3381

Zeocin™ ORF (no Met): bases 3382-3753

SV40 polyadenylation sequence: bases 3883-4013

pUC origin: bases 4396-5069 (opposite strand)

Ampicillin resistance gene: bases 5214-6074 (opposite strand)