



**Comments for pIZ/V5-His/CAT
3552 nucleotides**

- OpIE2 promoter: bases 4-552
- CAT ORF: bases 616-1272
- V5 epitope: bases 1339-1380
- 6xHis tag: bases 1390-1407
- OpIE2 Reverse priming site: bases 1417-1442
- OpIE2 polyadenylation sequence: bases 1425-1554
- pUC origin: bases 1623-2296
- OpIE2 promoter: bases 2341-2889
- EM7 promoter: bases 2907-2984
- Zeocin™ resistance marker (ORF): bases 2985-3359
- SV40 early polyadenylation sequence: bases 3423-3552