




Latex



Nitrile

## >> Features

- Available in powder free latex and powder free nitrile
- Textured fingertips for an improved grip  
Beaded cuff design helps users don gloves efficiently
- Extra length for increased protection  
Biocompatibility tested (ISO 10993)
- Non-sterile



Nitrile material promotes better puncture



Extra-long cuffs

Cat #	Material	Size	Color	Unit	Total
97-1111	Latex	XS	Natural Off-White	100 /Pack, 10 Packs/Case	1000 /Case
97-1112	Latex	S	Natural Off-White	100 /Pack, 10 Packs/Case	1000 /Case
97-1113	Latex	M	Natural Off-White	100 /Pack, 10 Packs/Case	1000 /Case
97-1114	Latex	L	Natural Off-White	100 /Pack, 10 Packs/Case	1000 /Case
97-1115	Latex	XL	Natural Off-White	100 /Pack, 10 Packs/Case	1000 /Case
97-6111	Nitrile	XS	Dark Blue	100 /Pack, 10 Packs/Case	1000 /Case
97-6112	Nitrile	S	Dark Blue	100 /Pack, 10 Packs/Case	1000 /Case
97-6113	Nitrile	M	Dark Blue	100 /Pack, 10 Packs/Case	1000 /Case
97-6114	Nitrile	L	Dark Blue	100 /Pack, 10 Packs/Case	1000 /Case
97-6115	Nitrile	XL	Dark Blue	100 /Pack, 10 Packs/Case	1000 /Case